



TECHNOPATH
CLINICAL DIAGNOSTICS

**SEROLOGY
CONTROLS**

Multichem® ID
INFECTIOUS DISEASE QC

Multichem® ID-COVID19Neg

Negative Third-Party Quality Control for
IgG and IgM Antibodies to SARS-CoV-2



Features

- Ready to use liquid stable controls
- Technical and scientific support

Specifications

- 36 month shelf life at 2°C to 8°C.
- 30 day open vial stability at 2°C to 8°C.

| Product | Description | Config. | Part Code |
|-------------------------|------------------|---------|-----------|
| Multichem ID-COVID19Neg | Negative Control | 4 x 4mL | CVN200N |
| Multichem ID-COVID19Neg | Negative Control | 4 x 2mL | CVN100N |

Analytes

Negative for SARS-CoV-2 IgG
Negative for SARS-CoV-2 IgM

For further information see Instructions For Use.

THE QUALITY CONTROL COMPANY

TECHNOPATH CLINICAL DIAGNOSTICS
www.technopathcd.com
multichemcovid@technopathcd.com

MM_IDCOVID19Neg_Digi_R2