



Date 25 January 2024

Technopath Clinical Diagnostics
Technopath Life Sciences Park
Fort Henry
Ballina
Co. Tipperary
V94 FF1P
IRELAND

Tel: +353 61 525 700

LGC Clinical Diagnostics (part of the LGC Group) are a global leader in the manufacture of the Multichem® range of Internal Quality Controls. Our products are used by thousands of laboratories in over 100 countries globally.

We manufacture our Multichem range of controls in our ISO 13485 accredited facility.

The term “third party” is used to describe a quality control product that helps provide an independent assessment of a diagnostic device or method.

Multichem controls are third-party, independent controls that are designed to be run on multiple systems from various instrument providers.

LGC Clinical Diagnostics has agreements in place with various instrument manufacturers, that allows them to distribute Multichem controls to their customers globally. The instrument manufacturer has no input into the production of Multichem controls, they are produced independently of the calibrator, reagent and instrument, delivering unbiased performance assessment. They are true third-party controls, meeting regulatory & accreditation requirements.

This independence enables the Multichem control material to closely mimic the performance of patient samples, and in doing so, provide an unbiased, independent assessment of analytical performance across platforms.

Sincerely,

Bernd Hass
Senior VP of Quality and Regulatory Affairs
LGC Clinical Diagnostics
Bernd.hass@lgcgroup.com